

manugraphy®

Accurate & reliable hand analysis

manugraphy® enables the analysis of the hand function at highest quality level.

Scan the **pressure distribution** to get a reliable and accurate **analysis of the hand.**


Hand pressure platform

manugraphy® key features:


- collect pressure and force data during static & dynamic grasping and holding
- work with reliable, individually calibrated, capacitive sensors to evaluate anatomical segments of the hand individually
- follow pre-defined protocols and create a report within seconds, automatically
- quickly integrate the sensor in your lab or medical environment and sync with other systems



Technical information



manugraphy® software features



novel GmbH (Global, GER)
Ganghoferstrasse 68b, 80339 Munich
tel: +49 (89) 417767-0
e-mail: info@novel.de
web: www.novel.de
copyright © novel GmbH - August 2026

novel electronics inc. (North America)
3367 Babcock Blvd, Suite 101
Pittsburgh, PA 15237
tel: +1 (412) 755-0200
e-mail: novelinc@novelusa.com
web: www.novelusa.com